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
## *e-corder*

*Models 210, 410, 821 & 1621*

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
# Glossary

This glossary covers terms used in this manual.

**AC coupling.** When AC coupling is chosen, a 0.1 Hz high-pass filter before the first amplification stage which can be used to remove signal drift.

**ADC** (analog—to—digital convertor). A device that converts analog information into some corresponding digital voltage or current.

**amplitude.** The maximum vertical distance of a periodic wave from the zero or mean position about which the wave oscillates.

**analog.** An analog signal varies continuously over time, rather than changing in discrete steps.

**analog input.** This refers to the connectors on the front of the e-corder marked **Input**. These inputs are designed to accept up to  $-10$  V. Inputs can be either single—ended or differential.

**analog output.** This refers to the connectors on the front of the e-corder marked **Output**. The analog output provides a software-controlled variable output ( $-10$  V) that can be used with applications either directly as a stimulator, or to control peripheral devices.

**analysis mode.** When the e-corder is not connected to the computer, then Chart or Scope software can be used to analyze and manipulate existing files if the analysis option is chosen.

**BNC** (bayonet nut connector). A sort of cable or connector; a BNC—to—BNC cable connects two BNC connectors.

**bridge transducer.** A type of transducer using a Wheatstone bridge circuit. In its basic form, the bridge consists of four two—terminal elements (usually strain gauges) connected to form a quadrilateral. An excitation source is connected across one diagonal, and the transducer output is taken across the other.

**bus.** A data—carrying electrical pathway.

**Chart.** The Chart software is supplied with the e-corder and emulates a multi—channel chart recorder.

**connector.** A plug, socket, jack, or port used to connect one electronic device to another (via a cable).

**CPU** (central processing unit). A hardware device that performs logical and arithmetical operations

on data as specified in the instructions: the heart of most computers.

**DAC** (digital—to—analog convertor). A device that converts digital information into some corresponding analog voltage or current.

**DC offset.** The amount of DC (direct current) voltage present at the output of an amplifier when zero voltage is applied to the input; or the amount of DC voltage present in a transducer in its equilibrium state.

**differential input.** Input using both positive and negative inputs on an e-corder. The recorded signal is the difference between the positive and negative input voltages: if both were fed exactly the same signal, zero would result. Can reduce the noise from long leads.

**digital.** Varying in discrete steps. A digital signal increments in discrete steps rather than varying continuously. The size of the steps is determined by the resolution (16 bit, 12 bit etc.) and the full scale range of the input channel.

**DIN** (Deutsche Industrie Norm). An electronic standard.

**eDAQ Amp.** Software controlled preamplifiers that can be used with your e-corder. eDAQ Amps are recognized automatically by the e-corder system and seamlessly integrated into its applications, operating under full software control.

**envelope.** The overall shape of a signal, outlined by the minimum and maximum recorded values. Often used to display quickly changing signals.

**excitation.** The voltage supplied to a transducer or other device in order to power it.

**external trigger.** The input connector on the front of the e-corder marked Trigger. This lets one start recording from an external source. The trigger level (the voltage needed to have an effect) depends on the hardware and cannot be changed. The /20 series e-corders can also be triggered by contact closure if this is set up in the software.

**filter.** A mechanism whereby the frequency components of the signal are altered. Low pass filters remove high frequency (noise) components. High pass filters remove low frequency (drift) components.

**frequency.** The number of cycles per to occur in a waveform per unit. Frequency is usually expressed in the SI unit of Hertz, Hz (cycles per second,  $s^{-1}$ ), but is also sometimes expressed in beats per minute (bpm,  $min^{-1}$ ), revolutions per minute (rpm,  $min^{-1}$ ), cycles per hour (cph,  $h^{-1}$ ), etc.

**frequency response.** The bandwidth in which a circuit passes a signal without significant attenuation. A low pass filter's frequency response is the frequency where the signal is attenuated by 3 decibels, that is the frequency in filtered signal has an amplitude of  $0.707 (1/\sqrt{2})$  of the frequency in the input signal.

**gain.** The amount of amplification of a signal.

**half—bridge transducer** A bridge transducer only using half of the full-bridge circuit. It consists of two elements of equal value with an excitation voltage applied across them. The output of the transducer is taken at the junction of the two elements.

**Hertz (Hz).** See frequency.

**high pass filter (HPF).** See filter.

**I<sup>2</sup>C.** The I<sup>2</sup>C connection is used by the e-corder to control eDAQ Amps.

**IEC.** International Electrotechnical Commission.

**low pass filter (LPF).** See filter .

**PCI** (peripheral component interconnect). An architecture for connecting peripheral devices (such as USB cards) to computers.

**Pod connector.** A special 8-pin DIN connector on some e-corder inputs providing differential or single-sided connection, and they can also provide –5 V power.

**Pods.** Miniature preamplifiers that connect to the e-corder Pod connectors.

**port.** A socket where you plug in a cable for connection to a network or a peripheral device.

**range.** The full scale input range is region over which signals can be recorded. The e-corder has ranges of –2 mV to –10 V, in 12 steps. (Range is inversely proportional to gain, the extent of amplification.)

**Scope.** Scope software is supplied with the e-corder. It emulates a two—channel storage oscilloscope or XYT plotter.

**serial.** A connection protocol for sending information sequentially, one bit at a time, over a single wire.

**transducer.** A physical device that converts one energy form into another (usually an electrical signal). For example transducers can provide an electrical signal related to the amount of measured force, displacement, temperature, pressure, pH, etc. Ideally the signal varies in a proportional

manner to the measured parameter, but this is not always so.


**trigger.** A signal such as a voltage pulse, used to determine when sampling will begin. Sampling can be made to begin when the trigger level is reached, after it, or even prior to it. See also external trigger.

**TTL** (transistor—transistor logic). A family of integrated circuits with bipolar circuit logic, used in computers and related devices. TTL is also a standard for interconnecting such ICs, defining the voltages used to represent logical zeroes and ones (binary 0 and 1).

**USB.** Universal serial bus.

**waveform.** The shape of a wave; a graph of a wave s amplitude versus time.

# Index

## A

AC coupling 33  
accuracy 23  
ADC 21, 33  
analog input 11, 21  
analog output 13, 22  
anti—aliasing filter 10  
auxiliary ground 16

## B

back panel 15  
bandwidth 10  
BNC 33  
bridge 33  
Bridge Amp 4  
bus 33

## C

Chart 3  
computer requirements  
    Macintosh 3  
    Windows 3  
Conductivity Pod 4  
connection to the computer 19  
connector 34  
CPU 34

## D

DAC 22, 34  
DC drift 10  
differential input 34  
digital-to-analog converter 22  
DIN 34  
DIN connector 12

## E

e-corder  
    specifications 29  
    system 2  
eDAQ Amps  
    Bridge Amp 4  
    GP Amp 4  
    pH Amp 4  
    Picostat 4  
    Potentiostat 4  
    QuadStat 4  
eDAQ Pods  
    Conductivity Pod 4  
    pH Pod 4  
    RTD Pod 4  
    Thermocouple Pod 4  
eDAQ web site 5  
envelope 34  
excitation 34  
external trigger 14, 22

## F

filter 34  
front panel 8

## G

gain 34  
GP Amp 4  
ground connection 16

## H

half—bridge 35

## I

I2C port 15  
input amplifier 10

## N

New Hardware wizard 28

## O

offset 34

## P

PCI 35  
pH Amp 4  
pH Pod 4  
Picostat 4  
Pod connector 12, 35  
Pods 4  
Potentiostat 4  
Power indicator 9, 19  
power requirements 18  
problems and solutions 27

## Q

QuadStat 4

## R

RS485 serial port 16  
RTD Pod 4

## S


Scope 3  
self-test 18  
serial port 16  
signal accuracy 23  
specifications 29  
Status indicator 9, 19

## T

Thermocouple Pod 4  
transducer 35  
trigger 14, 34, 35  
Trigger indicator 9  
TTL 35

## U

USB 35  
cables 20  
connection 19  
hubs 20  
port 16  
user modification 21

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Copyright © October 2004